

ITK Test Harness: Guidance for the Test Message Pack

Value	Display	Guidance and Conformance	Required	Test ID/ Comments
10001	Handling Specification Error	<p>A generic error code which gives a minimum level of assurance that systems can share the minimum information relating to Handling Specification faults.</p> <p>The code '10001' is mandatory to be supported and returned to sender. Note a more specific code must be returned where one exists for the error.</p>	MUST	ITK-TOC-TEST-10001 This test message will return the 10001 response code.
10002	Infrastructure Level Response Value - Processing Error	<p>The handling specification flag for an infrastructure level response is present but cannot be processed. For example, it may be unreadable or contain an incorrect value. This code is mandatory so must be supported.</p> <p>If the InfAckRequested value received is not 'true' or 'false' the response code '10002' must be returned to sender.</p>	MUST	ITK-TOC-TEST-10002-1 ITK-TOC-TEST-10002-2 These test messages will return the 10002 response code.
10003	Business Level Response Value - Processing Error	<p>The handling specification flag for the business level response is present but cannot be processed. For example, it may be unreadable or contain an incorrect value. This code is mandatory so must be supported.</p> <p>If the BusAckRequested value received is not 'true' or 'false' the response code '10003' must be returned to sender.</p>	MUST	ITK-TOC-TEST-10003 This test message will return the 10003 response code.
10004	Message Definition Value – Processing Error	<p>The handling specification value for MessageDefinition is present but cannot be processed. For example, may be unreadable or contain an incorrect value. This may also be returned when the message type is not supported (known) by the receiving system. This code must be supported and the following values are applicable:</p>	MUST	ITK-TOC-TEST-10004 This test message will return the 10004 response code.
		<p>Transfer of Care – Emergency Care eDischarge Implementation Guide 2.6.0-beta.</p> <p>The value must be: 'https://fhir.nhs.uk/STU3/MessageDefinition/ITK - EC-eDischarge-MessageDefinition-4'</p>		
		<p>Transfer of Care – Inpatient (Acute) Discharge eDischarge Implementation Guide 2.6.0-beta.</p> <p>The value must be: 'https://fhir.nhs.uk/STU3/MessageDefinition/ITK - eDischarge-MessageDefinition-4'</p>		
		<p>Transfer of Care – Mental Health eDischarge Implementation Guide 2.6.0-beta.</p> <p>The value must be:</p>		

ITK Test Harness: Guidance for the Test Message Pack

		'https://fhir.nhs.uk/STU3/MessageDefinition/ITK-MH-eDischarge-MessageDefinition-4'		
		Transfer of Care – Outpatient Letter Implementation Guide 2.6.0-beta. The value must be: 'https://fhir.nhs.uk/STU3/MessageDefinition/ITK-OPL-MessageDefinition-4'		
		Digital Medication - DM Immunization Implementation Guide 1.2.0-alpha 'https://fhir.nhs.uk/STU3/MessageDefinition/ITK-DM-Immunization-MessageDefinition-1'		
		Digital Medication - DM Immunization Implementation Guide 1.2.0-alpha 'https://fhir.nhs.uk/STU3/MessageDefinition/ITK-DM-Immunization-MessageDefinition-1' If any other value is received the response code '10004' must be returned unless is an incorrect version, then follow the guidance below for response code '10005' [Note the GP writeback message string is TBA]		
10005	Message Definition Version Value – Processing Error	The handling specification for Message Definition is present but the version is not supported by the receiving system. If the string received is correct but the version is incorrect, the response code '10005' must be returned to sender	MUST	ITK-TOC-TEST-10005 The response is triggered by a dummy value in the message. e.g. sending Practitioner.name.given.value="10005" This test message will return the 10005 response code.
10007	Sender Reference Value - Processing Error	The handling specification string for Sender Reference is present but cannot be processed. For example, may be unreadable, contain an incorrect value or the use of Sender Reference is not supported by receiving system. For the current use cases if any value other than the SenderReference default value of "None" is received the response code '10007' must be returned to sender.		ITK-TOC-TEST-10007 The response is triggered by a dummy value in the message. e.g. sending Practitioner.name.given.value="10007" This test message will return the 10007 response code.
10008	Handling Specification Business Rule Error	The Handling Specification usage does not match business rules for included Payload. For example, an acknowledgement flag defined as mandatory by the Payload specification is missing.	Not applicable for current use cases, may be	

ITK Test Harness: Guidance for the Test Message Pack

		This error code should not be used for current use cases as there are no defined process for the correct identification of these rules.	required for future use cases	
10009	Unreadable message received	A message has been received that is either corrupted or malformed and cannot be read by the receiving system. The response code '10009' must be supported and returned to sender.	MUST	ITK-TOC-TEST-10009 This test message will return the 10009 response code.
10010	Recipient Type - Processing Error	The handling specification for Recipient Type is present but cannot be processed. For example, may be unreadable or contain an incorrect value. The response code '10010' must be supported and returned to sender.	MUST	No test message
20013	Success	The Message has been processed successfully. A response will be returned stating the fact. However, the message may result in a further business level response, if the business level response request flag has been sent to "true". The response code '20013' must be supported and must be returned to sender.	MUST	ITK-TOC-TEST-20013 The message uses test data for the patient's NHS number as a trigger. This test message will return the 20013 response code.
20001	Unrecognised Recipient Person	The person referred to as the recipient in the ITK3 MessageHeader (when present) is not recognised. The response code '20001' must be supported and must be returned to sender.	MUST	ITK-TOC-TEST-20001 The message uses test GP data to trigger the response. This test message will return the 20001 response code.
20002	Unrecognised Recipient Organisation	The organization referred to as the recipient in the ITK3 MessageHeader (when present) is not recognised. The response code '20002' must be supported and returned to sender.	MUST	ITK-TOC-TEST-20002 The message uses test GP Practice data to trigger the response. e.g. receiving Organization ODS identifier.value This test message will return the 20002 response code.
20003	Unrecognised Sender	The organization or person referred to as the sender in the ITK3 MessageHeader is not recognised. Note: This code should not be used where the domain makes use of the "GP look-up" functionality in MESH. This error code should not be used for current use cases due to the potential use of GP look-up.	Not applicable for current use cases	
20004	Non Approved file	The Receiving system has received an attached file whose file type is not approved for the business	Not applicable	

ITK Test Harness: Guidance for the Test Message Pack

	type received as an attachment	domain. This error code should not be used for current use cases due to the fact no business domain has defined their non-approved attachment types.	for current use cases, may be required in future use cases	
20005	Unsupported file type received as an attachment	<p>The Receiving system has received an attached file which it does not support.</p> <p>Note: all receiving systems must support at least .pdf file extensions.</p> <p>The response code '20005' must be supported and returned to sender.</p>	MUST	<p>ITK-TOC-TEST-20005</p> <p>The response is triggered by a dummy value in the message.</p> <p>e.g. sending Practitioner.name.given.value="20005"</p> <p>This test message will return the 20005 response code.</p>
20006	ITK3 header validation failure	<p>The ITK3 header resources or elements are not correct or understandable. The ITK3 Bundle or ITK3 MessageHeader fail conformance against the FHIR schemas:</p> <p>http://hl7.org/fhir/STU3/fhir-all-xsd.zip</p> <p>or HAPI validation using the approved version of the HAPI libraries*</p> <p>*This will be the version used by the ITK3 Test Harness.</p>	MUST	<p>ITK-TOC-TEST-20006</p> <p>This test message will return the 20006 response code.</p>

ITK Test Harness: Guidance for the Test Message Pack

20007	Duplicate Message received	<p>Bundle with this message identifier has already been processed. A Payload with this message identifier has already been received and processed by this recipient. A message with this message header.id has already been received and processed.</p> <p>As stated in the FHIR standard:- An incoming message contains two identifiers: the Bundle.id and the MessageHeader.id. Each time a new message is created, it SHALL be assigned an identifier (MessageHeader.id) that is unique within that message stream. Note that since message streams are often merged with other streams, it is recommended that the identifier should be globally unique. This can be achieved by using a UUID or an OID. Each time a message is sent, the Bundle.id should be changed to a new value.</p> <p>See http://hl7.org/fhir/STU3/messaging.html</p> <p>The response code '20007' must be supported and returned to sender, and the message deemed to be a duplicate must be rejected regardless of the content of the payload.</p>	MUST	<p>ITK-TOC-TEST-20007</p> <p>The response is triggered by a dummy value in the message.</p> <p>e.g. sending Practitioner.name.given.value="20007"</p> <p>This test message will return the 20007 response code.</p>
20008	Duplicate Document received	<p>A Bundle with this document identifier has already been processed. A Payload with this document identifier has already been received and processed by this recipient.</p> <p>As stated in the FHIR standard:-</p> <p>The document identifier ((is) mandatory). This is found in <i>Bundle.identifier</i> and is globally unique for this instance of the document, and is never re-used, including for other documents derived from the same composition.</p> <p>See http://hl7.org/fhir/STU3/documents.html</p> <p>The response code '20008' must be supported and returned to sender, and the document deemed to be a duplicate must be rejected regardless of the content of the payload.</p>	MUST	<p>ITK-TOC-TEST-20008</p> <p>The response is triggered by a dummy value in the message.</p> <p>e.g. sending Practitioner.name.given.value="20008"</p> <p>This test message will return the 20008 response code.</p>

ITK Test Harness: Guidance for the Test Message Pack

20009	Payload validation failure	<p>The payload is not correct or understandable. The payload fails conformance test against the FHIR schemas:</p> <p>http://hl7.org/fhir/STU3/fhir-all-xsd.zip</p> <p>or HAPI validation using the approved version of the HAPI libraries*</p> <p>*This will be the version used by the ITK3 Test Harness.</p> <p>The response code '20009' must be supported and returned to sender.</p>	MUST	<p>ITK-TOC-TEST-20009</p> <p>The response is triggered by a dummy value in the message.</p> <p>e.g. sending Practitioner.name.given.value="20009"</p> <p>This test message will return the 20009 response code.</p>
20010	Unrecognised Payload Recipient Organisation	<p>The Recipient Organisation identified in the Payload, is not supported by this End Point (Receiving System).</p> <p>The response code '20010' must be supported and returned to sender.</p>	MUST	No test message
20011	Unrecognised Payload Recipient Person	<p>The Recipient person identified in the Payload, is not supported by this End Point (Receiving System).</p> <p>The response code '20011' must be supported and returned to sender.</p>	MUST	No test message
20012	Unauthorised Sender	<p>The Receiving system identified in the Payload is configured to reject messages from unauthorised senders. This code should not be used where the domain makes use of the "GP look-up" functionality in MESH.</p> <p>The response code '20012' should not be used</p>	N/A	
20014	Unable to Process Replacement Document	<p>A replacement document was received, but the receiver could not process the new document correctly and has therefore marked the new and original documents as "bad" on its system.</p> <p>The response code '20014' must be supported and returned to sender.</p>	MUST	No test message

ITK Test Harness: Guidance for the Test Message Pack

30001	Patient known here. (e.g. Patient is registered here)	The response code '30001' must be supported and returned to sender.	MUST	<p>ITK-TOC-TEST-30001</p> <p>The business response is triggered by test data (Patient's NHS number) carried in Patient.identifier.value.</p> <p>This test message will return the 30001 response code.</p>
30002	Patient not known here. (aka 'patient record not present in system')	The response code '30002' must be supported and returned to sender.	MUST	<p>ITK-TOC-TEST-30002</p> <p>The business response is triggered by test data (Patient's NHS number) carried in Patient.identifier.value.</p> <p>This test message will return the 30002 response code.</p>
30003	Patient no longer at this clinical setting	The response code '30003' must be supported and returned to sender.	MUST	<p>ITK-TOC-TEST-30003</p> <p>The business response is triggered by test data (Patient's NHS number) carried in Patient.identifier.value.</p> <p>This test message will return the 30003 response code.</p>